

BT2000 PROGRAMMABLE TIMER

Reproducible Timer For Burkard Flow Injection Valves



- ❑ Low Cost
- ❑ Simple operation
- ❑ Easily programmable
- ❑ Precise control

Introduction

The BT2000 is a programme protected versatile timer with hand operated pencil switch for loop filling flow injection. Suitable for Selective Ion and Flow Injection analysis the timer is an ideal option for low sample numbers when no interfaced automatic sample changer is required.

Timing ranges include: 000.1secs to 999.9secs and 00.01secs to 99.99secs. The timer has a panel mounted start button in addition to the external switch and is fitted with internal lithium batteries for initial programming. The unit is supplied with mains and interface cable.

Ordering Information:

Product Code	Description
BS02634	BT2000 Programmable Timer

The company reserves the right to make changes to the specification without notice.

